

Product Datasheet

Human Immunodeficiency Virus 1 (Group O gp41) Antigen, Recombinant V9001

Product Name	Human Immunodeficiency Virus 1 (Group O gp41) Antigen, Recombinant
Catalog Number	V9001
Species/Host	E. coli
Tested Applications	WB, ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	20 mM Na-carbonate, pH 9.6
Specificity	Group O variant, gp41
Molecular Weight	26.7 kDa