

## Product Datasheet

### Human Immunodeficiency Virus 1 (subtype C p24) Antigen, Recombinant V9002

<b>Product Name</b>	Human Immunodeficiency Virus 1 (subtype C p24) Antigen, Recombinant
<b>Catalog Number</b>	V9002
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	LF, ELISA, WB
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	His tag
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	PBS, pH 7.4
<b>Specificity</b>	Subtype C clone, p24
<b>Molecular Weight</b>	17.0 kDa