

## Product Datasheet

### Treponema pallidum (TnpA) Antigen, Recombinant V9004

<b>Product Name</b>	Treponema pallidum (TnpA) Antigen, Recombinant
<b>Catalog Number</b>	V9004
<b>Species/Host</b>	E. coli
<b>Tested Applications</b>	WB, ELISA, LF
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified recombinant antigen
<b>Storage</b>	≤ -20°C
<b>Note</b>	His tag
<b>Purity</b>	~95% by SDS-PAGE
<b>Buffer</b>	20 mM Na-carbonate, pH 10.0
<b>Specificity</b>	TnpA
<b>Molecular Weight</b>	16.4 kDa