

Product Datasheet

Treponema pallidum (p17) Antigen, Recombinant V9005

Product Name	Treponema pallidum (p17) Antigen, Recombinant
Catalog Number	V9005
Species/Host	E. coli
Tested Applications	WB, LF, ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	20mM Carbonate
Specificity	TpN17
Molecular Weight	16.4 kDa