

Product Datasheet

Hepatitis C virus (core) Antigen, Recombinant V9007

Product Name	Hepatitis C virus (core) Antigen, Recombinant
Catalog Number	V9007
Species/Host	E. coli
Tested Applications	WB, ELISA, LF
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	His tag
Purity	~95% by SDS-PAGE
Buffer	Phosphate buffered saline pH 7.4
Specificity	Core protein
Molecular Weight	48.7 kDa