

Product Datasheet

Oropouche Virus (OROV-NP) Antibody HM1405

Product Name	Oropouche Virus (OROV-NP) Antibody
Catalog Number	HM1405
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Oropouche virus
Alternative Names	OPV, OPOV
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% NaCl
Specificity	Oropouche virus nucleoprotein (NP)
Matched Pairs	HM1407, HM1406 or HM1405 (detection)
Preservatives	0.05% Sodium Azide