

Product Datasheet

Oropouche Virus (OROV-NP) Antibody HM1407

Product Name	Oropouche Virus (OROV-NP) Antibody
Catalog Number	HM1407
Description	Polyclonal Antibody, IgG Fraction
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Oropouche virus
Alternative Names	OPV, OPOV
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% NaCl
Specificity	Oropouche virus nucleoprotein (NP)
Matched Pairs	HM1406, HM1405 or HM1403 (detection)
Preservatives	0.05% Sodium Azide