

Product Datasheet

Alpha-1-Antitrypsin, (A1AT) Human Plasma 106-11-1

Product Name	Alpha-1-Antitrypsin, (A1AT) Human Plasma
Tested Applications	Control, Calibration
Brand	Lee Biosolutions
Concentration	Reported in mg per vial
Storage	−20 °C
Note	Reconstitute in DI water
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Buffer	Lyophilized from sodium phosphate buffer
Molecular Weight	~54,000 Da
Shelf Life	3 years