

## Product Datasheet

### MDMA-BSA Antigen 170056

<b>Product Name</b>	MDMA-BSA Antigen
<b>Catalog Number</b>	170056
<b>Description</b>	Synthetic MDMA, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Store at -20 °C. Short term storage ( $\leq$ 2 weeks): +2–8 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with a Catalog code P82-99-51A. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled Water
<b>Matched Pairs</b>	100906 (detection)
<b>Shelf Life</b>	60 months