



Date: 2017-12-14

PRODUCT SPECIFICATIONS

Name Anti-h Insulin 9350

Specificity Antibody recognizes human insulin

Description Monoclonal antibody raised against human Insulin. *In vitro* cultured mouse monoclonal

antibody.

Product code 700013

Product buffer solution 0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN₃ as a preservative

Shelf life and storage 60 months from manufacturing at 2-8 °C

Analyte description Insulin is a peptide hormone produced in the pancreas by the beta cells. It is formed as

a precursor molecule preproinsulin which is cleaved into proinsulin and further to mature insulin. Insulin is important for carbohydrate and fat metabolism and its

deficiency leads to diabetes mellitus.

PARAMETERS TESTED ON EACH LOT

Product appearance Liquid

Product concentration > 1.0 mg/ml

Immunoreactivity Passed in ELISA when Anti-h Insulin antibody 9350 was used as a detection antibody

IEF Profile 6.8 – 7.2

Purity ≥ 95 %

PARAMETERS DETERMINED DURING PRODUCT DEVELOPMENT

Subclass IgG₁

Association rate constant N/D

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





Cross-reactivities N/D

Epitope N/D

Pair recommendationsCAPTURE ANTIBODYDETECTION ANTIBODY

9351 9350

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore

these recommendations are only indicative.

Product stability TEMPERATURE, TIME RESULT

N/D

Miscellaneous

References -