



2021-02-19

Product specifications

Name Anti-h PTH (Parathyroid hormone) 10750

Specificity Antibody recognizes human parathyroid hormone (PTH)

Description Monoclonal antibody raised against human Parathyroid Hormone (PTH). In vitro cultured

mouse monoclonal antibody.

Product code 700024

Product buffer solution 0.1~M~PBS, pH 7.4, 0.9~%~NaCl, $0.05~\%~NaN_3$ as a preservative

Shelf life and storage 60 months from manufacturing at 2-8 °C

Subclass IgG_1

Analyte description Parathyroid hormone is a peptide hormone secreted by the parathyroid glands that regulates

the blood calcium and phosphate levels. PTH is routinely monitored for people with chronic kidney disease or who are on dialysis. Measurement of intact PTH is used to diagnose primary

hyperparathyroidism and hypercalcaemia.

Parameters tested on each lot

Product appearance Liquid

Product concentration > 1.0 mg/ml

Passed in RIA when Anti-h PTH antibody 10750 was used as a detection antibody **Immunoreactivity**

IEF Profile Not Determined (N/D)

N/D Purity

Kinetic parameters

Association rate constant N/D

Dissociation rate constant N/D

Affinity constant N/D

Determination method

Determination antigen





2021-02-19

N/D Cross-reactivities

Epitope Amino acids 57-61

Pair recommendations

Platforms tested RIA

Antigens tested N/D

Product stability TEMPERATURE, TIME **RESULT**

N/D

Miscellaneous

References