



Product specifications

Name Anti-h ApoA1 2002 SPTN-5

Specificity Antibody recognizes human apolipoprotein A1 (ApoA1)

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components

Product code 100264

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 36 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Apolipoprotein A1 (ApoA1) is a major component of the high-density lipoprotein complex

(HDL or "good cholesterol"). ApoA1 helps to clear cholesterol from arteries. Decreased serum HDL cholesterol levels have been reported to correlate with increased risk of coronary artery disease. However, ApoA1 has been suggested as a better discrimination of coronary artery

disease than HDL.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.7–7.7

Purity ≥ 95 %

Kinetic parameters

Association rate constant 8 x 10⁵ 1/Ms

Dissociation rate constant 8 x 10⁻⁴ 1/s

Affinity constant $K_A = 1 \times 10^9 \text{ 1/M}; K_D = 1 \times 10^{-9} \text{ M}$

Determination method SPR analysis (ProteOn XPR36)

Determination antigen Human Apolipoprotein A1, Chemicon (Cat ALP10, Lot 0507005923)





2021-03-31

Cross-reactivities Not Determined (N/D)

Epitope N/D

Pair recommendations

		DETECTION	
		2001	2002
CAPTURE	2001	-	+
	2002	+	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA

Antigens tested Recombinant ApoA1, Medix Biochemica, 710012, 710033, 710032.

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days Precipitation observed -20 °C, 21 days Precipitation observed

+4 °C, 21 days OK +25 °C, 21 days OK +35 °C, 21 days OK +45 °C, 3 days OK

+45 °C, 7 days Reduced immunoreactivity

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -