



## **Product specifications**

Name Anti-h CRP 6402 SPTN-5

Specificity Antibody recognizes human C-reactive protein (CRP), binding is Ca<sup>2+</sup> independent

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components.

Product code 100145

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 24 months from manufacturing at 2–8 °C

Subclass IgG<sub>1</sub>

Analyte description CRP is a member of the class of acute-phase reactants, as its levels rise dramatically during

inflammatory processes occurring in the body. CRP rises up to 50,000-fold in acute inflammation, such as infection. It rises above normal limits within 6 hours, and peaks at 48 hours. Its half-life is constant, and therefore its level is mainly determined by the rate of production (and hence the severity of the precipitating cause). CRP is used mainly as a marker of inflammation. Measuring and charting CRP values can prove useful in determining disease

progress or the effectiveness of treatments.

## Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.6–7.7

Purity ≥ 95 %

## Kinetic parameters

Association rate constant 7.3 x 10<sup>5</sup> 1/Ms

Dissociation rate constant 1.2 x 10<sup>-4</sup> 1/s

Affinity constant  $K_A = 6.2 \times 10^9 \text{ 1/M}$ ;  $K_D = 1.6 \times 10^{-10} \text{ M} (= 0.16 \text{ nM})$ 

Determination method SPR analysis (ProteOn XPR36)

Determination antigen Human CRP, Scripps (Cat C0124, Lot 2200701)





Cross-reactivities

Does not recognize serum amyloid-P complex

**Epitope** 

Not Determined (N/D)

Pair recommendations

		DETECTION				
		6402	6403	6404	6405	6407
CAPTURE	6402	+	+	+	+	+
	6403	+	+	+	+	-
	6404	+	+	+	+	+
	6405	+	-	+	+	+
	6407	+	-	+	+	+

Following pairs are especially recommended for the below mentioned assays:

LF: 6402 (membrane) – 6402 (particles), 6402 – 6404, 6402 – 6405 IT: 6405 – 6405, hs-CRP (high sensitivity): 6404 – 6407, 6405 – 6407

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested

FIA, LF, IT

Antigens tested

Native CRP antigen, Lee Biosolutions 140-11A and 140-12, recombinant CRP antigen, Lee Biosolutions 140-11R.

**Product stability** 

TEMPERATURE, TIME	RESULT
-70 °C, 21 days	OK
-20 °C, 21 days	OK
+4 °C, 21 days	OK
+25 °C, 21 days	OK
+35 °C, 21 days	OK
+45 °C, 7 days	OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous

Note that CRP is a pentameric protein. The same antibody can thus be used both as a capture and a detection antibody.

References

Käpyaho, K., Welin, M.-G., Tanner, P., Kärkkäinen, T. and Weber, T. (1989) Rapid determination of C-reactive protein by enzyme immunoassay using two monoclonal antibodies. Scand. J. Clin. Lab. Invest. 49:389-393