



Product specifications

Name Anti-h MPO 1702 SPTN-5

Specificity Antibody recognizes human myeloperoxidase

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components

Product code 100267

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 36 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Myeloperoxidase (MPO) is a peroxidase enzyme (EC 1.11.1.7) most abundantly present in

neutrophil granulocytes. The 150-kDa MPO is a dimer consisting of two 15-kDa light chains and two variable-weight glycosylated heavy chains bound to a prosthetic heme group. Myeloperoxidase (MPO) is an enzymatic mediator of several inflammatory cascades and higher serum levels have been associated with increased risk of adverse cardiovascular events. As a result, MPO and its downstream inflammatory pathways represent attractive targets for both prognostic and therapeutic intervention in the prophylaxis of atherosclerotic

cardiovascular disease.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.7–7.5

Purity ≥ 95 %

Kinetic parameters

Association rate constant 8 x 10⁵ 1/Ms

Dissociation rate constant 2 x 10⁻³ 1/s

Affinity constant $K_D = 2 \times 10^{-9} \text{ M}$; $K_A = 5 \times 10^8 \text{ 1/M}$

Determination method SPR analysis (ProteOn XPR36)

Determination antigen Myeloperoxidase, Lee Biosolutions (Cat 426-10, Lot L44340)





2021-03-31

Cross-reactivities Not Determined (N/D)

Epitope N/D

Pair recommendations

		DETECTION		
		1701	1702	1703
CAPTURE	1701	-	+	+
	1702	-	-	+
	1703	+	+	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA

Antigens tested Native MPO antigen Lee Biosolutions 426-10.

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days OK
-20 °C, 21 days OK
+4 °C, 21 days OK
+25 °C, 21 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -