



Product specifications

Name Anti-h PCT 4003 SPTN-5

Specificity Antibody recognizes human procalcitonin (PCT)

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components

Product code 100562

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage Unspecified, storage at 2–8 °C

Subclass IgG₁

Analyte description Procalcitonin (PCT) is a precursor protein of the hormone calcitonin, produced by the thyroid

gland. It is a 116 amino acid protein that is cleaved intracellularly to form three peptides: N-terminal procalcitonin, calcitonin, and katacalcin. In healthy individuals, only the calcitonin polypeptide is secreted into the bloodstream, and PCT serum levels are very low. However, in patients with systemic infections, and especially sepsis, serum PCT levels rise rapidly.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.4–6.7

Purity ≥ 95 %

Kinetic parameters

Association rate constant 4.8 x 10⁴ 1/Ms

Dissociation rate constant 2.2 x 10⁻⁴ 1/s

Affinity constant $K_A = 2.2. \times 10^8 \text{ 1/M}; K_D = 4.6 \times 10^{-9} \text{ M} (=4.6 \text{ nM})$

Determination method SPR analysis (ProteOn XPR36)

Determination antigen PCT, Medix Biochemica (Cat 610080)





2021-03-31

Cross-reactivities

Does not recognize canine procalcitonin.

Epitope

Anti-h PCT 4003 antibody recognizes the katacalcin domain of procalcitonin, i.e. the region Asp(96) - Asn(116), as calculated from the peptide sequence of the mature PCT protein, without the signal sequence.

Pair recommendations

		DETECTION						
		4003	4004	4005	4006	4008	4051	
CAPTURE	4003	-	+	+	-	-	+	
	4004	+	ı	-	+	+	-	
	4005	+	ī	-	+	+	-	
	4006	-	+	+	-	-	+	
	4008	-	+	+	-	-	+	
	4051	+	-	-	+	+	-	

Following pairs are especially recommended for the below mentioned assays: LF: 4006 (membrane) – 4004 (particles) and 4006 (membrane) – 4051 (particles)

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested

FIA, LF, CLIA

Antigens tested

Recombinant PCT antigen, Medix Biochemica 610080.

Product stability

TEMPERATURE, TIME	RESULT
-70 °C, 21 days	OK
-20 °C, 21 days	OK
+4 °C, 21 days	OK
+35 °C, 21 days	OK
+45 °C, 7 days	OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous

-

References

-