



## **Product specifications**

Name	Anti-Rotavirus VP6 11101 SPTN-5
Specificity	Antibody recognizes Rotavirus VP6
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components
Product code	100715
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_{3}$ as a preservative
Shelf life and storage	18 months from manufacturing at 2–8 °C
Subclass	IgG <sub>1</sub>
Analyte description	Rotavirus infection is the leading cause of severe diarrhea and dehydration in young children. Rotavirus viral protein 6 (VP6), which is highly conserved and the most abundant protein on the viral capsid, is a good target for diagnostic use in Rotavirus immunoassays.

## Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80–120 % compared to the reference sample in an FIA test
IEF Profile	6.8–7.5
Purity	≥ 95 %

## **Kinetic parameters**

Association rate constant	Not Determined (N/D)
Dissociation rate constant	N/D
Affinity constant	N/D
Determination method	-
Determination antigen	-

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. Medix/MAB and Medix/Antigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.





Cross-reactivities	N/D		
Epitope	N/D		
Pair recommendations	CAPTURE ANTIBODY 11101	DETECTION ANTIBODY 11101	
	Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.		
Platforms tested	FIA		
Antigens tested	Recombinant Rotavirus VP6 antigen, Medix Biochemica 610145.		
Product stability		RESULT OK OK OK OK OK OK outfer to see whether different temperatures ition of the antibody. Please note that the shelf me stability testing at 2–8 °C in the product	
Miscellaneous	-		
References	-		

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI14631532

Legal disclaimer Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.