



2022-01-12

Product specifications

Name Anti-h ST2 10205 SPTN-5

Specificity Antibody recognizes human soluble ST2

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components.

Product code 100684

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 36 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description ST2 is an interleukin-1 family receptor expressed in the heart. ST2 protein has two isoforms: a

soluble (sST2) and membrane-bound form. Patients with heart failure and elevated ST2

protein levels in the blood are at risk for heart failure progression.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 5.6–6.3

Purity ≥ 95 %

Kinetic parameters

Association rate constant 2.8 x 10⁵ 1/Ms

Dissociation rate constant 4.8 x 10⁻⁴ 1/s

Affinity constant $K_A = 5.8 \times 10^8 \text{ 1/M}$; $K_D = 2.0 \times 10^{-9} \text{ M}$ (= 2.0 nM)

Determination method BLI (Octet RED96e)

Determination antigen Human ST2, Acro Biosystems IL1-H5229





2022-01-12

Cross-reactivities

Epitope Recognizes free sST2

N/D

Pair recommendations

		DETECTION						
		10201	10202	10203	10204	10205	10206	10207
CAPTURE	10201	-	-	-	+	+	+	+
	10202	-	-	-	+	+	+	+
	10203	-	-	-	+	+	+	+
	10204	+	+	+	ı	+	+	-
	10205	+	+	+	+	-	-	-
	10206	+	+	+	+	-	-	-
	10207	+	+	+	-	-	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA, CLIA

Antigens tested Recombinant ST2 antigen, Medix Biochemica, 710020

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days OK
-20 °C, 21 days OK
+4 °C, 21 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -