



## **Product specifications**

Name Anti-Thyroxine 6902 SPTN-5

Specificity Antibody recognizes human thyroxine

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from animal-derived

components.

Product code 100801

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 18 months from manufacturing at 2–8 °C

Subclass IgG<sub>1</sub>

Analyte description Thyroxine, or 3,5,3',5'-tetraiodothyronine (often abbreviated as T4) is the major hormone

secreted by the thyroid gland. T4 is transported in blood, with 99.95 % of the secreted T4 being protein bound, principally to thyroxine-binding globulin (TBG). T4 is involved in controlling the rate of metabolic processes in the body and influencing physical development. Thyroxine is a prohormone and a reservoir for the active thyroid hormone triiodothyronine

(T3) which is about four times more potent.

## Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

 ${\color{red}Immunoreactivity} {\color{gray}80-120\,\%\,compared\,to\,the\,reference\,sample\,in\,an\,FIA\,test}$ 

IEF Profile 6.7–7.6

Purity ≥ 95 %

## Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





Cross-reactivities T3 (3,3',5-triiodothyronine) 1.5 % (Sigma, Cat T2752) rT3 (3,3',5'-triiodothyronine, reverse T3) 1.5 % (Sigma, Cat T0281)

3,5-diiodothyronine < 0.1 % (Sigma, Cat D0629)

Epitope N/D

Pair recommendations CAPTURE ANTIBODY DETECTION ANTIBODY

Platforms tested FIA

Antigens tested Thyroxine (T4) antigen, BSA conjugate, Lee Biosolutions 581-10

Product stability TEMPERATURE, TIME RESULT

 -70 °C, 21 days
 OK

 -20 °C, 21 days
 OK

 +4 °C, 21 days
 OK

 +35 °C, 21 days
 OK

 +45 °C, 7 days
 OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -