

Product specifications

Name	Ketamine-BSA Antigen
Description	Synthetic Ketamine, BSA conjugated
Amino acid sequence	-
Product host	-
Product code	170037
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Ketamine is a dissociative drug and anesthetic with hallucinogenic properties.
Product concentration	Lot-specific
Miscellaneous	This product was previously available with a Catalog code P82-99-31A. It has been produced using the same conjugation protocol.
Matched pair	100902, P02-99-31M-P, HM577 or PP2111F (detection)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.

