



Product specifications

Name Propoxyphene (PPX)-BSA Antigen

Description Synthetic Propoxyphene, BSA conjugated

Amino acid sequence -

Product host -

Product code 170032

Product formulation Liquid; Antigen-BSA conjugate

Product buffer solution Distilled water

Reconstitution -

Shelf life and storage Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid

repeated freezing and thawing.

Analyte description Propoxyphene is an opioid used to treat pain.

Product concentration Lot-specific

Miscellaneous This product was previously available with a Catalog codes P001A and CB501. It has been

produced using the same conjugation protocol.

Matched pair 100895, P001M or CB302 (detection)