



Date: 2018-03-19

PRODUCT SPECIFICATIONS

Name Recombinant Calprotectin S100A9, 1000 μg

Description Recombinant human S100 calcium binding protein A9 antigen is produced by *E.coli*

expression system and the target gene encoding Thr2-Pro114 is expressed. Predicted

molecular weight: 13 kDa.

Amino acid sequence TCKMSQLERNIETIINTFHQYSVKLGHPDTLNQGEFKELVRKDLQNFLKKENKNEKVIEHIMEDL

DTNADKQLSFEEFIMLMARLTWASHEKMHEGDEGPGHHHKPGLGEGTP

Production host Escherichia coli (E. coli)

Product code 710036

Product formulation Lyophilized

Product buffer solution PBS, pH 7.4

Reconstitution Reconstitute lyophilized protein with 2000 µl of deionized water

Shelf life and storage -20 °C


Analyte description Calprotectin belongs to the family of \$100 proteins and is composed of two subunits

S100A8 and S100A9. Increased plasma concentrations are found in response to infections and inflammation. Faecal calprotectin measurement is used for detection of gastrointestinal inflammation or neoplasia and for assessment of inflammatory bowel

disease activity and response to treatment.

Product concentration 0.5 mg/ml when reconstituted with 2000 μl of deionized water

Purity SDS-PAGE > 90 %



MK R

N: Non-reducing R: Reducing

MK: Protein Molecular Weight Marker

Reactivity with MedixMABs Anti-h Calprotectin 3403: +

Anti-h Calprotectin 3404: -

Anti-h Calprotectin 3405: -

Anti-h Calprotectin 3406: -

Anti-h Calprotectin 3407: -