



## **Product specifications**

Name Recombinant Calprotectin antigen, 100 μg

Description Recombinant human calprotectin heterotetramer (S100A8/S100A9) protein with an N-

terminal histidine tag in the S100A8 subunit. Predicted molecular weights: S100A8 12 kDa,

S100A9 13 kDa.

Amino acid sequence Not Determined (N/D)

Product host Escherichia coli (E. coli)

Product code 610061

Product formulation Lyophilized

Product buffer solution Lyophilized from 20 mM Tris-HCl, 100 mM NaCl, 1 mM EDTA, pH 8.0; containing 6 % sucrose as

a stabilizer

Reconstitution Reconstitute the lyophilized protein with deionized water

Shelf life and storage 24 months for lyophilized product, storage at 2-8 °C. After reconstitution 1 month at 2-8 °C

and 1 month at -20 °C. Reconstituted product is recommended to be stored in plastic (PP)

containers.


Analyte description Calprotectin belongs to the family of S100 proteins and is composed of two subunits S100A8

and S100A9. Increased plasma concentrations are found in response to infections and inflammation. Fecal calprotectin measurement is used for detection of gastrointestinal inflammation or neoplasia and for assessment of inflammatory bowel disease activity and

response to treatment.

Product concentration N/A

Purity Capillary electrophoresis (CE-SDS)



Reactivity with MedixMABs Anti-h Calprotectin 3403: +

Anti-h Calprotectin 3404: +

Anti-h Calprotectin 3405: +

Anti-h Calprotectin 3406: +

Anti-h Calprotectin 3407: +

**Oy Medix Biochemica Ab** Klovinpellontie 3 FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg.no. FI146315320

