



Date: 2020-11-27

Product specifications

Name Recombinant IL-6, 500 μg

Description Recombinant human interleukin-6 antigen is produced by *E.coli* expression system and the

target gene encoding Pro29-Met212 is expressed. Predicted molecular weight: 19 kDa.

Amino acid sequence PVPPGEDSKDVAAPHRQPLTSSERIDKQIRYILDGISALRKETCNKSNMCESSKEALAENNLNLP

KMAEKDGCFQSGFNEETCLVKIITGLLEFEVYLEYLQNRFESSEEQARAVQMSTKVLIQFLQKKA

KNLDAITTPDPTTNASLLTKLQAQNQWLQDMTTHLILRSFKEFLQSSLRALRQM

Product host Escherichia coli (E. coli)

Product code 710039

Product formulation Lyophilized

Product buffer solution 20 mM Hac-NaAc, pH 5.0, 150 mM NaCl

Reconstitution Reconstitute the lyophilized protein with 1000 µl of deionized water

Shelf life and storage Unspecified, storage at -20 °C

Analyte description IL-6 is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. It is

secreted by T cells and macrophages to stimulate immune response to trauma, especially burns or other tissue damage leading to inflammation. IL-6 is relevant to many disease processes such as diabetes, atherosclerosis, depression, Alzheimer's Disease, systemic lupus

erythematosus, cancer, and rheumatoid arthritis.

Product concentration 0.5 mg/ml when reconstituted with 1000 µl of deionized water

Purity SDS-PAGE > 90 %

Reacitivity with MedixMABs Anti-h IL-6 2703: +

Anti-h IL-6 2704: + Anti-h IL-6 2706: + Anti-h IL-6 2707: +