



Date: 2018-08-22

PRODUCT SPECIFICATIONS

Name Recombinant ST2, 1000 μg

Description Recombinant human ST2, also known as interleukin-1 receptor-like 1, antigen is

produced by mammalian expression system and the target gene encoding Lys19-Phe328 is expressed with a 6His tag at the C-terminus. Predicted molecular weight: 60

kDa.

Amino acid sequence KFSKQSWGLENEALIVRCPRQGKPSYTVDWYYSQTNKSIPTQERNRVFASGQLLKFLPAAVADSG

IYTCIVRSPTFNRTGYANVTIYKKQSDCNVPDYLMYSTVSGSEKNSKIYCPTIDLYNWTAPLEWF KNCQALQGSRYRAHKSFLVIDNVMTEDAGDYTCKFIHNENGANYSVTATRSFTVKDEQGFSLFPV IGAPAQNEIKEVEIGKNANLTCSACFGKGTQFLAAVLWQLNGTKITDFGEPRIQQEEGQNQSFSN GLACLDMVLRIADVKEEDLLLQYDCLALNLHGLRRHTVRLSRKNPSKECFVDHHHHHH

Production host Human cells

Product code 710046

Product formulation Lyophilized

Product buffer 20 mM Na-phosphate, pH 7.4, 150 mM NaCl

Reconstitution Reconstitute lyophilized protein with 2000 µl of deionized water


Shelf life and storage -20 °C

Analyte description ST-2 is an interleukin-1 family receptor expressed in the heart. ST2 protein has two isoforms: a soluble (sST2) and membrane- bound form. Patients with heart failure and

elevated sST2 protein levels in the blood are at risk for heart failure progression.

Product concentration 0.5 mg/ml when reconstituted with 2000 µl of deionized water

Purity SDS-PAGE > 90 %



N: Non-reducing R: Reducing

MK: Protein Molecular Weight Marker

Reactivity with MedixMAB

Anti-h ST2 10201: +
Anti-h ST2 10202: +
Anti-h ST2 10203: +
Anti-h ST2 10204: +
Anti-h ST2 10205: +
Anti-h ST2 10206: +
Anti-h ST2 10207: +