

Ketamine-BSA Antigen

| | |
|--------------------------------|---|
| Catalog number | 170037 |
| Description | Synthetic Ketamine, BSA conjugated |
| Amino acid sequence | - |
| Product host | - |
| Product formulation | Liquid; Antigen-BSA conjugate |
| Product buffer solution | Distilled water |
| Reconstitution | - |
| Shelf life and storage | Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing. |
| Analyte description | Ketamine is a dissociative and anesthetic drug with hallucinogenic properties. |
| Product concentration | Lot-specific |
| Matched pairs | 100902, P02-99-31M-P, HM577 or PP2111F (detection) |
| Miscellaneous | This product was previously available with a Catalog code P82-99-31A. It has been produced using the same conjugation protocol. |



Legal disclaimer