

Methylphenidate-BSA Antigen

Catalog number	170039
Description	Synthetic Methylphenidate, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Methylphenidate is a central nervous system stimulant drug used to treat Attention Deficient Hyperactivity Disorder (ADHD) and narcolepsy. Abuse risk increases when administered in excessive doses and via intravenous and intranasal routes.
Product concentration	Lot-specific
Matched pairs	P03-99-01R-IF, P03-99-02R-IF, P03-99-03R-IF (detection)
Miscellaneous	This product was previously available with a Catalog code P83-99-01A. It has been produced using the same conjugation protocol.

**Legal disclaimer**