

Propoxyphene (PPX)-BSA Antigen

Catalog number	170032
Description	Synthetic Propoxyphene, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Propoxyphene is an opioid used to treat pain.
Product concentration	Lot-specific
Matched pairs	100895, P001M or CB302 (detection)
Miscellaneous	This product was previously available with a Catalog codes P001A and CB501. It has been produced using the same conjugation protocol.



Legal disclaimer