


Recombinant *H. pylori* OMP, 100 µg

Catalog number	610140
Description	Recombinant <i>Helicobacter pylori</i> outer membrane protein (OMP) with a histidine tag. Predicted molecular weight 60 kDa.
Amino acid sequence	Not Determined (N/D)
Product host	<i>Escherichia coli</i> (<i>E. coli</i>)
Product formulation	Liquid, slightly yellow
Product buffer solution	50 mM Sodium phosphate pH 7.0, 0.095 % NaN ₃ as a preservative
Reconstitution	-
Shelf life and storage	Unspecified, storage at 2-8 °C, do not freeze
Analyte description	<i>Helicobacter pylori</i> is a Gram-negative bacterium commonly found in the stomach. <i>H. pylori</i> infection increases the risk of developing peptic ulcers and gastric cancers and can be diagnosed with serological blood tests, breath tests and stool antigen tests.
Product concentration	Lot specific
Purity	Capillary electrophoresis (CE-SDS)



Antibodies tested Anti-*H. pylori* 7105: +

