

Tricyclic Antidepressant (TCA)-BSA Antigen

Catalog number	170030
Description	Synthetic Tricyclic Antidepressant, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	Shelf life unspecified. Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Analyte description	Tricyclic antidepressants are a class of medications used primarily for the treatment of depression and other mood disorders.
Product concentration	Lot-specific
Matched pairs	T002M or CB300 (detection).
Miscellaneous	This product was previously available with Catalog codes T001A and CB500. It has been produced using the same conjugation protocol.

**Legal disclaimer**