

SPECIFICATION

Aprotinin Bovine Lung

GENERAL INFORMATION

Catalogue Number	125-09	Storage Temperature	2-8°C
Lot Number	Reported	Storage Note	Avoid freeze-thaws
CAS Number	9087-70-1	Shipping Note	Cold packs
Molecular Weight	6511.44	Received By	_____
		Receipt Date	_____

PRODUCT INFORMATION

Source	Bovine Lung
Form	0.2 um filtered solution in 0.9% sodium chloride and 0.9% benzyl alcohol.
Appearance	Clear, colorless solution
Activity	6 - 16 TIU/mL
Protein	1-3mg protein/mL (Biuret)
Specific Activity	≥ 5 TIU/mg protein
Recertification	2 years from date of manufacture

HANDLING INFORMATION

Handling	Handle as a potentially hazardous substance.
----------	--

This product is intended for research, laboratory use, or further manufacturing. It is not for human use.

