Medix Biochemica

6-Monoacetylmorphine (6-MAM)-BSA Antigen

Catalog number	170031
Description	Synthetic 6-Monoacetylmorphine, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
Analyte description	Heroin (diacetylmorphine) is metabolized to 6-monoacetylmorphine and morphine. 6-monoacetylmorphine is a unique metabolite of heroin and its presence is definitive evidence of recent heroin use. Morphine is also a metabolite of several legal drugs, including codeine, and therefore is not a conclusive indication of heroin use.
Product concentration	Lot-specific
Matched pairs	JP091, JP092, JP093, HM445, HM446 or HM461 (detection).
Miscellaneous	This product was previously available with Catalog codes P83-99-41A and LA263. It has been produced using the same conjugation protocol.



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.