


NGAL Antigen, Recombinant

Catalog number	610012
Description	Recombinant human neutrophil gelatinase-associated lipocalin (NGAL) protein with a C-terminal histidine tag. Predicted molecular weight: 22 kDa.
Product host	<i>Escherichia coli (E. coli)</i>
Product formulation	Lyophilized
Product buffer solution	Lyophilized from Phosphate buffered saline, pH 7.4; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers
Reconstitution	Reconstitute the lyophilized protein with deionized water
Recertification	30 months
Storage	Storage at 2-8 °C
Product concentration	N/A
Purity	Capillary electrophoresis (CE-SDS)



Antibodies tested	Anti-h NGAL 4202: + Anti-h NGAL 4203: + Anti-h NGAL 4204: + Anti-h NGAL 4205: +
--------------------------	--



Legal disclaimer