

2023-04-20

Anti-h PTH (Parathyroid hormone) 10750

Product overview

Catalog number 700024

Specificity Antibody recognizes human parathyroid hormone (PTH)

Description Monoclonal antibody raised against human Parathyroid Hormone (PTH).

In vitro cultured mouse monoclonal antibody.

Product buffer solution 0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN₃ as a preservative

Shelf life and storage 60 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Parathyroid hormone is a peptide hormone secreted by the parathyroid

glands that regulates the blood calcium and phosphate levels. PTH is routinely monitored for people with chronic kidney disease or who are on

dialysis. Measurement of intact PTH is used to diagnose primary

hyperparathyroidism and hypercalcemia.

Parameters tested on each lot

Product appearance Liquid

Product concentration > 1.0 mg/ml

Immunoreactivity Passed in RIA when Anti-h PTH antibody 10750 was used as a detection

antibody

IEF Profile Not Determined (N/D)

Purity N/D

Kinetic parameters

Association rate constant N/D

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





Product specification ANTIBODY

2023-04-20

Cross-reactivities N/D

Epitope Amino acids 57-61

Pair recommendations

Platforms tested RIA

Antigens tested N/D

Product stability TEMPERATURE, TIME RESULT

N/D

Miscellaneous

References

