

## Product specification ANTIBODY

2023-05-24

## Anti-h CA15-3 R4406 SPTN-5

#### **Product overview**

Catalog number 140033

**Specificity** Antibody recognizes sialylated carbohydrate epitope on the MUC-1 protein

of cancer antigen 15-3

**Description** Recombinant mouse antibody, cultured *in vitro* under conditions free from

animal-derived components.

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN<sub>3</sub> as a preservative

Shelf life and storage Unspecified, storage at 2–8 °C

Subclass IgG<sub>1</sub>

**Analyte description** Cancer Antigen 15-3 is defined by the use of two monoclonal antibodies,

one specific for MUC-1 protein core and another antibody specific for a sialylated carbohydrate epitope on the MUC-1 protein. CA15-3 is a widely used serum marker in breast cancer patients for monitoring their response

to cancer therapy.

#### Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

**Product concentration** 5.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.6)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.7–7.2

**Purity** ≥ 95 %

### Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





# Product specification ANTIBODY

2023-05-24

Cross-reactivities

**Epitope** Antibody binds to sialylated carbohydrate epitope on the MUC-1 protein.

Pair recommendations

		DETECTION	
		4404	R4406
CAPTURE	4404	-	+
	R4406	+	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA

Antigens tested Native CA15-3 antigen Medix Biochemica 151-53A.

Product stability TEMPERATURE, TIME RESULT

N/D

-70 °C, 21 days OK -20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK

+35 °C, 21 days Charge alterations +45 °C, 7 days Charge alterations

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

**Miscellaneous** Antibody has similar epitope to 115D8.

References -

