

Product specification ANTIBODY

2022-11-17

Anti-h IgE 8516 SPTN-5

Product overview

Catalog number 100858

Specificity Antibody recognizes human immunoglobulin E

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from

animal-derived components.

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 24 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Immunoglobulin E (IgE) is a class of antibody that has only been found in

mammals. It plays an important role in allergy and defense against parasites. Although IgE is typically the least abundant isotype - blood serum IgE levels in a normal individual are only 0.05 % of the IgG concentration - it is capable of triggering the most powerful immune

reactions.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.1–7.1

Purity ≥ 95 %

Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant $K_A = 9.2 \times 10^9 \text{ 1/M}$

Determination method SPR analysis (Biacore)

Determination antigen IgE, Biodesign (Cat A101644, Lot 912737)





Product specification ANTIBODY

2022-11-17

Cross-reactivities Human IgA <0.32 %

Human IgG <0.32 % Human IgM <0.32 %

Epitope N/D

Pair recommendations

		DETECTION		
		8510	8516	8518
CAPTURE	8510	-	+	+
	8516	-	-	-
	8518	-	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA

Antigens tested N/D

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days N/D
-20 °C, 21 days OK
+4 °C, 21 days OK
+30 °C, 7 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -

