Medix Biochemica

Anti-h PG I 8016 SPRN-5

Product overview						
Catalog number	100868					
Specificity	Antibody recognizes human pepsinogen I					
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.					
Product buffer solution	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 $\%$ NaCl, 0.095 $\%$ NaN $_3$ as a preservative					
Shelf life and storage	36 months from manufacturing at 2–8 °C					
Subclass	lgG₁					
Analyte description	Pepsinogen is the pro-form of pepsin and is produced in the stomach by chief cells. The major part of pepsinogen is secreted into the gastric lumen but a small amount can be found in the blood. Alterations in the serum pepsinogen concentrations has been found with <i>Helicobacter pylori</i> (<i>H. pylori</i>) infections, peptic ulcer disease, gastritis, and gastric cancer. More precise analysis may be achieved by measuring the pepsinogen I/II ratio.					
Parameters tested on each lot						
Product appearance	Liquid, may turn slightly opaque during storage					
Product concentration	5.0 mg/ml (+/- 10 %)					
Immunoreactivity	80–120 % compared to the reference sample in an FIA test					
IEF Profile	5.2-8.2					
Purity	≥ 95 %					
Kinetic parameters						
Association rate constant	Not Determined (N/D)					
Dissociation rate constant	N/D					
Affinity constant	3 x 10 ¹⁰ 1/M					
Determination method	Radioimmunoassay (RIA)					
Determination antigen	PG I, Medix Biochemica in-house (TK 05)					



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Cross-reactivities	Human Pepsinogen II Porcine Pepsinogen Gastrin I Somatostatin Peptide YY Bombesin Cathepsin D			0.54 % (Medix Biochemica in-house, A69/96) 65.2 % (Sigma, Cat P-4656, Lot 82H8065) < 0.01 % (Sigma, Cat G-9020, Lot 110H48131) < 0.01 % (Sigma, Cat S-9129, Lot 44H03692) < 0.01 % (Sigma, Cat P-1306, Lot 73H49502) < 0.01 % (Sigma, Cat B-4272, Lot 95H58651) < 0.06 % (Sigma, Cat C-8696, Lot 46H0795)					
Epitope	N/D						1		
Pair recommendations			DETECTION						
			8003	8009	8015	8016	_		
	Ш	8003	-	+	+	+	-		
	TUF	8009	+	-	+	-	-		
	CAPTURE	8015	+	+	-	+	-		
	0	8016	+	+	-	-			
	Following pairs are especially recommended for the below mentioned assays: FIA: 8003 (capture) – 8009 (detection), 8003 – 8015, 8015 – 8003, 8016 – 8003 and 8016 – 8009 LF: 8003 (membrane) – 8009 (particles), 8003 – 8015, 8009 – 8015 and 8015 – 8009 IT: 8003 – 8016 Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.								
Platforms tested	FIA, LF, IT								
Antigens tested	Recombinant Pepsinogen I antigen, Medix Biochemica 610000 and native Pepsinogen I antigen Lee Biosolutions 440-50.								
Product stability	TEMPERATURE, TIMERESULT-70 °C, 21 daysN/D-20 °C, 21 daysOK+4 °C, 21 daysOK+30 °C, 21 daysOK+35 °C, 21 daysOK+45 °C, 7 daysOKStability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.						ion of the		
Miscellaneous	-								
References	-								



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.