

Product specification ANTIBODY

2023-02-10

Anti-h PRL 5602 SP-5

Product overview

Catalog number 100360

Specificity Antibody recognizes human prolactin (luteotropic hormone, or luteotropin)

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from

animal-derived components.

Product buffer solution 0.9 % NaCl, 0.095 % NaN₃ as a preservative

Shelf life and storage 24 months from manufacturing at 2–8 °C

Subclass IgG₁

Analyte description Prolactin (PRL) or luteotropic hormone (LTH) is a peptide hormone

primarily associated with lactation. Prolactin is usually measured when checking for pituitary tumors and the cause of: breast milk production that is not related to childbirth (galactorrhea), impotence, infertility, irregular

menstrual periods.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/-10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 5.7–6.7

Purity $\geq 95\%$

Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant 9.5 x 10⁹ 1/M

Determination method Radioimmunoassay (RIA)

Determination antigen PRL, UCB –Bioproducts (Cat i117/001, Lot 003)





Product specification ANTIBODY

2023-02-10

Cross-reactivities GH < 0.5 % (Boehringer, Cat 252964, Lot 1222303/Jul 84)

PL < 0.05 % (UCB-Bioproducts, Cat i020, Lot 007)

Epitope N/D

Pair recommendations

		DETECTION	
		5601	5602
CAPTURE	5601	-	+
	5602	+	-

Following pairs are especially recommended for the below mentioned assays:

LF: 5602 (membrane) – 5601 (particles)

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA, CLIA, LF

Antigens tested Native PRL antigen, Lee Biosolutions 996-41.

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days N/D
-20 °C, 21 days N/D
+4 °C, 21 days N/D
+35 °C, 7 days N/D
+35 °C, 21 days N/D
+45 °C, 3 days N/D
+45 °C, 7 days N/D

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References

Oy Medix Biochemica Ab

FI-02180 Espoo, Finland

Klovinpellontie 3,

