

Product specification ANTIBODY

2023-01-20

Anti-h SCC 11803 SPTN-5

Product overview

Catalog number 100878

Specificity Antibody recognizes human squamous cell carcinoma (SCC) antigen

Description Monoclonal mouse antibody, cultured *in vitro* under conditions free from

animal-derived components.

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9% NaCl, 0.095% NaN₃ as a preservative

Shelf life and storage Unspecified, storage at 2–8 °C

Subclass IgG₁

Analyte description Squamous cell carcinoma (SCC) antigen is a biomarker of squamous cell

carcinomas found in the lung, uterine cervix, esophagus, head, neck, anal canal, and skin. In clinical applications SCC antigen levels are measured from serum or plasma to monitor recurrence of the tumor and response to the therapy. This marker is not recommended for screening, since SCC antigen levels may also increase in non-malignant diseases including

severe inflammatory diseases, psoriasis, and renal failure.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 5.0 mg/ml (+/- 10 %)

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.1–6.5

Purity $\geq 95 \%$

Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





Product specification ANTIBODY

2023-01-20

Cross-reactivities N/D

Epitope N/D

Pair recommendations

		DETECTION			
		11801	11802	11803	11804
CAPTURE	11801	-	-	+	+
	11802	-	-	+	+
	11803	+	+	-	+
	11804	+	+	+	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA, CLIA

Antigens tested N/D

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days OK
-20 °C, 21 days OK
+4 °C, 21 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

Miscellaneous -

References -

