

## Cannabinoid (THC)-BSA Antigen

---

<b>Catalog number</b>	170044
<b>Description</b>	Synthetic Cannabinoid, BSA conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-BSA conjugate
<b>Product buffer solution</b>	Distilled water
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	60 months from manufacturing at -20 °C. Short term storage ( $\leq 2$ weeks): 2–8 °C. Avoid repeated freezing and thawing.
<b>Analyte description</b>	Tetrahydrocannabinol (THC) is a psychoactive cannabinoid.
<b>Product concentration</b>	Lot-specific
<b>Matched pairs</b>	100892, 100893, 100894, 100898 or P01-99-21M-P (detection)
<b>Miscellaneous</b>	This product was previously available with a Catalog code P81-99-21A. It has been produced using the same conjugation protocol.

**Legal disclaimer**