

## Product specification ANTIBODY

2024-01-17

# Cannabinoid (THC) Antibody

#### **Product overview**

Catalog number 100892

**Specificity** Antibody recognizes THC and its metabolites

**Description** Protein A purified monoclonal mouse antibody, cultured *in vitro* under

conditions free from animal-derived components.

Product buffer solution 0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN<sub>3</sub>

as a preservative

**Shelf life and storage** 18 months from manufacturing at 2–8 °C

Subclass IgG<sub>1</sub>

**Analyte description** Tetrahydrocannabinol (THC) is a psychoactive cannabinoid.

### Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

**Product concentration** ≥ 1.5 mg/ml

Immunoreactivity 80–120 % compared to the reference sample in an FIA test

IEF Profile 6.9–7.8

**Purity**  $\geq 95\%$ 

### Kinetic parameters

Association rate constant N/D

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -





## Product specification ANTIBODY

2024-01-17

**Cross-reactivities** Does not react with Opiates, Benzoylecgonine, Amphetamines,

Phencyclidine and Cannabidiol.

Epitope N/D

Pair recommendations 170044 (BSA conjugate as capture)

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and

therefore these recommendations are only indicative.

Platforms tested ELISA, LF

Antigens tested Please see Pair recommendations above

Product stability TEMPERATURE, TIME RESULT

N/D

Miscellaneous This product was previously available with Catalog codes P01-99-23M-P

and AG309. It has been produced using the same original mouse

monoclonal hybridoma cell line.

Immunogen: 8-THC-BSA

References -

