

2023-03-08

Cortisol-3-HRP Antigen

Catalog number 170053

Description Synthetic Cortisol, HRP conjugated

Amino acid sequence -

Product host -

Product formulation Liquid; Antigen-HRP conjugate

Product buffer solution 1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005%

thimerosal

Reconstitution -

Shelf life and storage Shelf life unspecified, storage at 2–8 °C. Do not freeze.

Analyte description Cortisol is a steroid hormone produced by the adrenal glands. It is involved

in many important functions of the body including regulating stress response, reducing inflammation, maintaining blood pressure and controlling blood sugar. Cortisol is often measured together with ACTH to diagnose disorders of the pituitary or adrenal glands such as Cushing's

syndrome or Addison's disease.

Product concentration -

Purity -

Antibodies tested 100904, P01-92-94M-P (capture)

Miscellaneous This product was previously available with catalog code P91-92-91H. It has

been produced using the same conjugation protocol.

