

## Cotinine-4-BSA Antigen

<b>Catalog number</b>	170047
<b>Description</b>	Synthetic Cotinine, BSA conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-BSA conjugate
<b>Product buffer solution</b>	Distilled Water
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	Shelf life unspecified. Short term storage ( $\leq 2$ weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Analyte description</b>	Cotinine is a biomarker of tobacco exposure and the main metabolite of nicotine.
<b>Product concentration</b>	Lot specific
<b>Purity</b>	-
<b>Matched pairs</b>	PP194IF, PP195IF, PP198IF (detection)
<b>Miscellaneous</b>	This product was previously available with catalog code P82-99-04A. It has been produced using the same conjugation protocol.

**Legal disclaimer**