

Cotinine-4-HRP Antigen

Catalog number	170049
Description	Synthetic Cotinine, HRP conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-HRP conjugate
Product buffer solution	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
Reconstitution	-
Shelf life and storage	Shelf life unspecified, storage at 2–8 °C. Do not freeze.
Analyte description	Cotinine is a biomarker of tobacco exposure and the main metabolite of nicotine.
Product concentration	-
Matched pairs	PP194IF, PP195IF, PP198IF (capture)
Miscellaneous	This product was previously available with a Catalog code P92-99-02H. It has been produced using the same conjugation protocol.



Legal disclaimer