

Fentanyl-BSA Antigen

Catalog number	170040
Description	Synthetic Fentanyl, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
Analyte description	Fentanyl is a synthetic opioid with 50–100 times the potency of morphine. It is used to treat severe pain. Fentanyl may be illicitly produced and mixed with other illegal drugs increasing the abuse and overdose risks.
Product concentration	Lot-specific
Matched pairs	P01-99-53R-IF, P01-99-54R-IF, P01-99-55R-IF (detection)
Miscellaneous	This product was previously available with a Catalog code P81-99-51A. It has been produced using the same conjugation protocol.

**Legal disclaimer**