

K2-BSA Antigen

Catalog number	170082
Description	Synthetic K2, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
Analyte description	K2 (also known as Spice) is a mixture of dried and shredded plant material sprayed with synthetic cannabinoid substances. It may also be found in liquid form. The effects of K2 are unpredictable and often dangerous.
Product concentration	Lot-specific
Matched pair	JP098, JP099 or JP100 (detection)
Miscellaneous	This product was previously available with a Catalog code P83-99-51A and LA288. It has been produced using the same conjugation protocol.

**Legal disclaimer**