

## Ketamine-BSA Antigen

---

|                                |   |
|--------------------------------|---|
| <b>Catalog number</b>          | 170060  |
| <b>Description</b>             | Synthetic Ketamine, BSA conjugated  |
| <b>Amino acid sequence</b>     | -   |
| <b>Product host</b>            | -   |
| <b>Product formulation</b>     | Liquid; Antigen-BSA conjugate   |
| <b>Product buffer solution</b> | Distilled water   |
| <b>Reconstitution</b>          | -   |
| <b>Shelf life and storage</b>  | 60 months from manufacturing at -20 °C. Short term storage ( $\leq 2$ weeks): 2–8 °C. Avoid repeated freezing and thawing.                      |
| <b>Analyte description</b>     | Ketamine is a dissociative and anesthetic drug with hallucinogenic properties. Ketamine's chemical structure is similar to phencyclidine (PCP). |
| <b>Product concentration</b>   | Lot-specific  |
| <b>Matched pairs</b>           | 100902 (detection)  |
| <b>Miscellaneous</b>           | This product was previously available with a Catalog code P82-99-32A. It has been produced using the same conjugation protocol.                 |

**Legal disclaimer**