

## Methaqualone-BSA Antigen

---

<b>Catalog number</b>	170072
<b>Description</b>	Synthetic Methaqualone, BSA conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-BSA conjugate
<b>Product buffer solution</b>	Distilled water
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
<b>Analyte description</b>	Methaqualone is a synthetic central nervous system depressant that acts as a sedative and hypnotic. It leads to drowsiness and lowered blood pressure. It is dangerous when used in excess or in combination with other depressants.
<b>Product concentration</b>	Lot-specific
<b>Matched pair</b>	P02-99-81R-IF, P02-99-82R-IF or P02-99-83R-IF (detection)
<b>Miscellaneous</b>	This product was previously available with a Catalog code P82-99-81A. It has been produced using the same conjugation protocol.

**Legal disclaimer**