

## Methylphenidate-BSA Antigen

---

<b>Catalog number</b>	170039
<b>Description</b>	Synthetic Methylphenidate, BSA conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-BSA conjugate
<b>Product buffer solution</b>	Distilled water
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	60 months from manufacturing at -20 °C. Short term storage ( $\leq 2$ weeks): 2–8 °C. Avoid repeated freezing and thawing.
<b>Analyte description</b>	Methylphenidate is a central nervous system stimulant drug used to treat Attention Deficient Hyperactivity Disorder (ADHD) and narcolepsy. Abuse risk increases when administered in excessive doses and via intravenous and intranasal routes.
<b>Product concentration</b>	Lot-specific
<b>Matched pairs</b>	P03-99-01R-IF, P03-99-02R-IF, P03-99-03R-IF (detection)
<b>Miscellaneous</b>	This product was previously available with a Catalog code P83-99-01A. It has been produced using the same conjugation protocol.

**Legal disclaimer**