

## Recombinant Calprotectin S100A8, 50 µg

Catalog number 710018

**Description** Recombinant human S100 calcium binding protein A8 antigen is produced

by E. coli expression system and the target gene encoding Met1-Glu93 is expressed with a 6His tag at the C-terminus. Predicted molecular weight:

14 kDa.

Amino acid sequence MLTELEKALNSIIDVYHKYSLIKGNFHAVYRDDLKKLLETECPQYIRKKGAD

VWFKELDINTDGA VNFQEFLILVIKMGVAAHKKSHEESHKEHHHHHH

**Product host** Escherichia coli (E. coli)

**Product formulation** Lyophilized

**Product buffer solution** 20 mM HEPES, pH 7.4, 150 mM NaCl

Reconstitution Reconstitute lyophilized protein with 100 µl of deionized water

-20 °C Storage at


**Analyte description** Calprotectin belongs to the family of S100 proteins and is composed of two

> subunits S100A8 and S100A9. Increased plasma concentrations are found in response to infections and inflammation. Fecal calprotectin measurement is used for detection of gastrointestinal inflammation or neoplasia and for assessment of inflammatory bowel disease activity and response to

treatment.

**Product concentration** 0.5 mg/ml when reconstituted with 100 µl of deionized water

**Purity** SDS-PAGE > 90 %



N: Non-reducing R: Reducing

MK: Protein Molecular Weight Marker

**Antibodies tested** 

Anti-h Calprotectin 3403: -Anti-h Calprotectin 3404: + Anti-h Calprotectin 3405: + Anti-h Calprotectin 3406: + Anti-h Calprotectin 3407: -

www.medixbiochemica.com

VAT reg. no. FI14631532

