

Recombinant Influenza A NP, 100 µg

Catalog number 610050

Description Recombinant influenza A nucleoprotein antigen (strain A/Hong Kong/1/1968

H3N2) with a histidine tag. Predicted molecular weight: 57 kDa.

Amino acid sequence Not Determined (N/D)

Product host Escherichia coli (E. coli)

Product formulation Liquid

Product buffer solution 50 mM Tris-HCl, 250 mM NaCl, 10 % glycerol, pH 8.0

Reconstitution -

Shelf life and storage Unspecified, storage at -20 °C


Analyte description Influenza A is a genus of the Orthomyxoviridae family of viruses. It includes

one viral species which causes influenza in birds and some mammals. Influenza A subtypes are labeled according to their hemagglutinin and neuraminidase protein variants. Influenza nucleoprotein (NP) binds to viral RNA, and these ribonucleoprotein complexes play a major role in viral transcription and replication. NP is highly conserved among influenza A virus strains, making NP a very suitable control antigen to be used in

influenza A immunological tests.

Product concentration Lot specific

Purity Capillary electrophoresis (CE-SDS)



Antibodies tested Anti-Influenza A 7304: +

Anti-Influenza A 7305: + Anti-Influenza A 7307: +

